Lyophilized Exosome Standards



Purified exosomes from various sources

Product Code: HBM-###-#####

Pioneering the EV Field

1- About Exosomes

Exosomes are EVs actively secreted by exocytosis by most living cells, typically with diameter reported between 40 and 120 nm. Exosome release occurs either constitutively or upon induction, under both normal and pathological conditions. Both quantity and molecular composition of released exosomes depend on the physiological state of the parental cells.

2- Lyophilized Exosomes Standards

Lyophilized Exosomes are purified by a combination of tangential flow filtration (TFF), size exclusion chromatography (SEC). Isolated vesicles are quantified and validated for total protein content, common marker expression (CD9, CD81,CD63), particle size distribution and concentration by NTA (Nanoparticles Tracking Analysis) with Zetaview analyzer (Particle Metrix). Lyophilization does not alter the stability of exosome proteins and nucleic acids, in comparison to other storage methods, including storage of fresh exosomes at -20°C. Lyophilized exosomes are easy to ship and stable for long term storage (up to 36 months).

3- Types of Exosome Standards available

- Lyophilized Exosome Standards from human Biofluids (plasma, serum, urine, saliva) of healthy donors.
- Lyophilized Exosomes from cell culture media (COLO1, BLCL21, HCT116, SK-N-SH, U87, PC3, BPH-1, DAUDI, A549, K562, mouse cell B16F10).
- Lyophilized Exosomes from Human Mesenchymal Stem Cells (MSC) from adipose tissue. Pool of 10 different donors.
- Lyophilized Exosome Standards size avalable: 100 μg and 30 μg, sold in packages of 2, 5 vials.

4- Procedure for Exosome Standards reconstitution


- Reconstitute Lyophilized Exosome Standard by adding deionized water, 100 µl for Lyophilized Standard 100 µg and 30 µl for Lyophilized Standard 30 µg, to get a final concentration of 1 µg/µL. Different volumes of deionized water for exosomes reconstitution can be choosen by the users in according with the desired final concentration. Resuspend exosomes pipetting the solution up and down 10-15 times, avoiding bubbles. Vortex the reconstituted standard for 60 seconds.
- Briefly centrifuge the tubes containing the standard to ensure that the solution is collected at the bottom of the tube. Pipette the solution up and down 10 times, avoiding the introduction of bubbles. After this step, the standard is ready to use.

5- Storage

- Lyophilized Exosomes can be stored for 36 months at 4°C.
- Reconstituted Exosome Standards are not suitable for long term conservation at room temperature; use them within 2 hours after reconstitution. The remaining reconstituted solution should be aliquoted into polypropylene vials (preferably low binding) and stored at -20°C for up to one month or at -80°C for up to six months. Strictly avoid repeated freeze-and-thaw cycles.

6- Performance

Lyophilization does not affect EV particle size distribution or biomarker expression compared to other storage methods (Fig 1). Exosomes stored for over 3 months at -20° C or over 1 year at -80° C showed a different size distribution profile, probably due to EV aggregation.



1. Western Blot comparison of exosomal markers on fresh (F), frozen (-20°C) and lyophilized exosomes (L). Particle size distribution of Exosomes stored lyophilized or frozen.

7- Application of Lyophilized Exosome Standards

7a- Marker detection by Western Blotting.


Reconstituted Exosomes can be directly lysed in Laemmli buffer, then loaded on the Electrophoresis gel.

Recommended quantity: 10-20 µg per line.

7b- Phenotyping by ELISA assay.

Reconstituted Exosomes can be loaded directly onto ELISA plate wells.

Recommended quantity: 10 - 20 µg per well.






7c- NTA in fluorescence and scattered mode.

Reconstituted Exosomes can be used for phenotyping assays by fluorescence NTA.

Recommended quantity: 10 µg of exosomes incubated with 1 µg* of antibody fluorophore conjugated or 2 µl* of membrane dye CMDR.




*Right exosome/dye ratio needs to be determined by the user.



7d- Phenotyping assays by FACS.


Reconstituted Exosomes can be used for profiling biomarkers by FACS analysis.

Recommended quantity: 5 µg of reconstituted Exosomes Standards for each test.



7e- Profiling of Exosome associated nucleic acids.

Reconstituded Exosomes can be used for profiling nucleic acids (RNAs and DNA) biomarkers. The lysis buffer can be added directly to the Exosome preparation.



Recommended quantity: 20-50 µg per reaction.

8- References

Islam, M. K. B., & Marcus, R. K. (2025). Effects of packing density and adsorption conditions on extracellular vesicle dynamic binding capacities for Capillary-Channeled polymer (C-CP) fiber columns. Journal of Chromatography A, 466068. https://doi.org/10.1016/j.chroma.2025.466068

Suresh, P. S., & Zhang, Q. (2025). Comprehensive Comparison of Methods for Isolation of Extracellular Vesicles from Human Plasma. Journal of Proteome Research. https://doi.org/10.1021/acs.jproteome.5c00149

Olsen, D. A., Hansen, T. F., Brandslund, I., & Madsen, J. S. (2025). Development of a high-sensitivity assay for direct quantitation of EGFr on extracellular vesicles in plasma. Journal of Immunological Methods, 113877. https://doi.org/10.1016/j.jim.2025.113877

Escudero-Cernuda, S., Eiro, N., Fraile, M., Vizoso, F. J., Fernández-Colomer, B., & Fernández-Sánchez, M. L. (2025). Limitations and challenges in the characterization of extracellular vesicles from stem cells and serum. Microchimica Acta, 192(5). https://doi.org/10.1007/s00604-025-07147-4

De Putter, M. P. P., Capuano, A., Numan, M., Hankemeier, T., & Shakalisava, Y. (2025). Enhanced separation of extracellular vesicles using capillary isotachophoresis with spacer compounds. Electrophoresis. https://doi.org/10.1002/elps.202400113

Zohouri, D., Mai, T. D., Reyre, M., Smadja, C., Krupova, Z., Talbot, L., & Taverna, M. (2025). Elucidation of extracellular vesicles behavior during capillary isoelectric focusing. Talanta, 128055. https://doi.org/10.1016/j.talanta.2025.128055

Kundu, S., Guo, J., Islam, Md. S., Rohokale, R., Jaiswal, M., & Guo, Z. (2025). A New Strategy to Functionalize Exosomes via Enzymatic Engineering of Surface Glycans and its Application to Profile Exosomal Glycans and Endocytosis. Advanced Science. https://doi.org/10.1002/advs.202415942

Zhou, X., Jaiswal, M., Shi, J., Guo, J., Kundu, S., Guo, Z., & Zeng, Y. (2024). Efficient enzymatic glycan engineering of extracellular vesicles using Nanomaterial-Interfaced microfluidics. ACS Applied Materials & Interfaces. https://doi.org/10.1021/acsami.4c20294

Druzhkova, T. A., Zhanina, M. Y., Vladimirova, E. E., Guekht, A. B., & Gulyaeva, N. V. (2024). Proteomic Spectrum of Serum Exosomes in Ischemic Stroke Patients Is Associated with Cognitive Impairment in the Post-Stroke Period. Biochemistry (Moscow), 89(9), 1595–1609. https://doi.org/10.1134/s0006297924090062

Escudero-Cernuda, S., Eiro, N., Fraile, M., Vizoso, F., Fernández-Colomer, B., & Sánchez, M. L. F. (n.d.). Characterization and quantification of human uterine cervical stem cells extracellular vesicles. https://papers.ssm.com/sol3/papers.cfm?abstract_id=4905805

Vu, V., Vu, C., Huang, C., Cheng, C., Pan, S., & Chen, W. (2024). Intermittent lysis on a single paper-based device to extract exosomal nucleic acid biomarkers from biological samples for downstream analysis. Microchimica Acta, 191(8). https://doi.org/10.1007/s00604-024-06566-z

Cho, Y., Cho, M. Y., Yoon, J., Hong, D. E., Lee, J., Park, H. S., Lee, H., Hong, K. S., Won-Kyu, L., Saehae, C., Song, S., & Noh, Y. (2024). Evaluation of unmodified human cell-derived extracellular vesicle mitochondrial deoxyribonucleic acid-based biodistribution in rodents. Journal of Extracellular Vesicles, 13(7). https://doi.org/10.1002/jev2.12489

Fujiwara, K., Takagi, Y., Tamura, M., Omura, M., Morimoto, K., Nakase, I., Tokonami, S., & Iida, T. (2023). Ultrafast sensitivity-controlled and specific detection of extracellular vesicles using optical force with antibody-modified microparticles in a microflow system. Nanoscale Horizons, 8(8), 1034–1042. https://doi.org/10.1039/d2nh00576j

Ng, C. T., Yip, W. K., Mohtarrudin, N., Jabar, M. F., & Seow, H. F. (2023). Increased tube formation and up-regulation of FGFR3 mRNA expression in microvascular endothelial cell by exosomes derived from \$W480-7. PubMed, 45(2), 195–204. https://pubmed.ncbi.nlm.nih.gov/37658529

Cheng, C., Hou, K., Hsu, C., & Chiang, L. (2023). Ultrasensitive and High-Resolution protein spatially decoding framework for tumor extracellular vesicles. Advanced Science, 11(3). https://doi.org/10.1002/advs.202304926

Jackson, K. K., Mata, C., & Marcus, R. K. (2023). A rapid capillary-channeled polymer (C-CP) fiber spin-down tip approach for the isolation of plant-derived extracellular vesicles (PDEVs) from 20 common fruit and vegetable sources. Talanta, 252, 123779.

Sandúa, A., Sanmamed, M. F., Rodríguez, M., Ancizu-Marckert, J., Gúrpide, A., Perez-Gracia, J. L., ... & González, Á. (2023). PSA reactivity in extracellular microvesicles to commercial immunoassays. Clinica Chimica Acta, 543, 117303. saliva

Vago, R., Radano, G., Zocco, D., & Zarovni, N. (2022). Urine stabilization and normalization strategies favor unbiased analysis of urinary EV content. Scientific reports, 12(1), 1-20.

Billotto, L. S., Jackson, K. K., & Marcus, R. K. (2022). Determination of the Loading Capacity and Recovery of Extracellular Vesicles Derived from Human Embryonic Kidney Cells and Urine Matrices on Capillary-Channeled Polymer (C-CP) Fiber Columns. Separations, 9(9), 251.

Otaki, M., Hirane, N., Natsume-Kitatani, Y., Nogami Itoh, M., Shindo, M., Kurebayashi, Y., & Nishimura, S. I. (2022). Mouse tissue glycome atlas 2022 highlights inter-organ variation in major N-glycan profiles. Scientific reports, 12(1), 1-15.

Kim, C., Han, B. H., Kim, D., Lee, G., Hong, C., Kang, J. Y., & Lee, K. (2022). Multi-target polydiacetylene liposome-based biosensor for improved exosome detection. Sensors and Actuators B: Chemical. 355, 131286.

Jackson, K. K., Powell, R. R., Bruce, T. F., & Marcus, R. K. (2022). Facile, generic capture and on-fiber differentiation of exosomes via confocal immunofluorescence microscopy using a capillary-channeled polymer fiber solid-phase extraction tip. Sensors & Diagnostics.

Riazanski, V., Mauleon, G., Lucas, K., Walker, S., Zimnicka, A. M., McGrath, J. L., & Nelson, D. J. (2022). Real time imaging of single extracellular vesicle pH regulation in a microfluidic cross-flow filtration platform. Communications biology, 5(1), 1-13.

Morasso, C., Ricciardi, A., Sproviero, D., Truffi, M., Albasini, S., Piccotti, F., ... & Corsi, F. (2022). Fast quantification of extracellular vesicles levels in early breast cancer patients by Single Molecule Detection Array (SiMoA). Breast cancer research and treatment, 192(1), 65-74.

Kruse, T., Schneider, S., Reger, L. N., Kampmann, M., & Reif, O. W. (2022). A novel approach for enumeration of extracellular vesicles from crude and purified cell culture samples. Engineering in Life Sciences.

Jackson, K. K., & Marcus, R. K. (2022). Rapid isolation and quantification of extracellular vesicles from suspension-adapted human embryonic kidney cells using capillary-channeled polymer fiber spin-down tips. Electrophoresis, 44(1–2), 190–202. https://doi.org/10.1002/elps.202200149

Choi, D. Y., Park, J. N., Paek, S. H., Choi, S. C., & Paek, S. H. (2022). Detecting early-stage malignant melanoma using a calcium switch-enriched exosome subpopulation containing tumor markers as a sample. Biosensors and Bioelectronics, 198, 113828.

Guo, K., Li, Z., Win, A., Coreas, R., Adkins, G. B., Cui, X., ... & Zhong, W. (2021). Calibration-free analysis of surface proteins on single extracellular vesicles enabled by DNA nanostructure. Biosensors and Bioelectronics, 192, 113502.

Li, S., Zhu, L., Zhu, L., Mei, X., & Xu, W. (2021). A sandwich-based evanescent wave fluorescent biosensor for simple, real-time exosome detection. Biosensors and Bioelectronics, 113902.

Miljkovic-Licina, M., Arraud, N., Zahra, A. D., Ropraz, P., & Matthes, T. (2022). Quantification and Phenotypic Characterization of Extracellular Vesicles from Patients with Acute Myeloid and B-Cell Lymphoblastic Leukemia. Cancers, 14(1),

El Ouahabi, O., Salim, H., Pero-Gascon, R., & Benavente, F. (2020). A simple method for the analysis of extracellular vesicles enriched for exosomes from human serum by capillary electrophoresis with ultraviolet diode array detection. Journal of Chromatography A, 461752.

Ridolfi, A., Brucale, M., Montis, C., Caselli, L., Paolini, L., Borup, A., ... & Nejsum, P. (2020). AFM-based high-throughput nanomechanical screening of single extracellular vesicles. Analytical chemistry, 92(15), 10274-10282.

Han, P., Lai, A., Salomon, C., & Ivanovski, S. (2020). Detection of Salivary Small Extracellular Vesicles Associated Inflammatory Cytokines Gene Methylation in Gingivitis. International journal of molecular sciences, 21(15), 5273.

